

AMIS[®] BIKINI

TAKE YOUR ANTERIOR APPROACH TO THE NEXT LEVEL

DEDICATED INSTRUMENTATION
TAILORED EDUCATION



Surgical Technique

Joint

Spine

Sports Med

ADDENDUM

TABLE OF CONTENTS

1. INTRODUCTION	4
1.1 The AMIS Bikini: the evolution of the AMIS approach	4
2. PATIENT SELECTION	4
3. PREOPERATIVE PLANNING	4
4. THE AMIS BIKINI APPROACH	5
4.1 Mobile Leg Positioner	5
4.2 Patient positioning	5
4.3 Surgical exposure	5
5. FEMORAL NECK OSTEOTOMY	6
6. ACETABULAR STAGE	6
7. FEMORAL STAGE	7
8. REDUCTION	8
9. WOUND CLOSURE	8
10. REFERENCES	9

1. INTRODUCTION

This document describes the surgical technique to perform a total or partial Hip arthroplasty using the bikini incision. This document represents an addendum to the AMIS surgical technique (ref. 99.98.12 or ref. 99.98.12US).


The AMIS Bikini is recommended for surgeons who have properly mastered the standard AMIS technique.

Dedicated AMIS Bikini instruments have been designed to enable the surgeon to perform an optimized and reproducible implantation via a bikini incision, further minimizing the risk of soft tissue damage respect to the traditional AMIS instruments.

NOTE: Please note that the instrument trays needed to perform a standard AMIS approach are always needed to perform the bikini incision technique in addition to the specific AMIS Bikini Add-on complement.

1.1 THE AMIS BIKINI: THE EVOLUTION OF THE AMIS APPROACH

The AMIS Bikini represents an evolution of the AMIS approach that follows the same intermuscular pathway with an inguinal skin incision performed within the skin fold of the "bikini line", or frontal groin crease. It results in an aesthetically pleasing cosmetic scar along the natural "Langer's lines" of the skin, that can be hidden when wearing underwear or swimwear (e.g., a bikini).



The AMIS Bikini combines all the benefits of the AMIS technique with the additional advantage of the groin crease incision which are:

- Better cosmetic/aesthetic appearance.^{1,2,3,4}
- Lower rates of delayed wound healing (also indicated patient with a high BMI).^{1,2,3,4}



2. PATIENT SELECTION

NOTE: For the patient selection not described in this addendum, please refer to the AMIS surgical technique (ref. 99.98.12 or ref. 99.98.12US).

The inguinal skin incision is particularly indicated in case of high BMI patient because of the lower rates of delayed wound healing.^{1,2,3,4}

3. PREOPERATIVE PLANNING

NOTE: For the preoperative planning please refer to the AMIS surgical technique (ref. 99.98.12 or ref. 99.98.12US).

4. THE AMIS BIKINI APPROACH

4.1 MOBILE LEG POSITIONER

NOTE: For the AMIS Mobile Leg Positioner please refer to the AMIS surgical technique (ref. 99.98.12 or ref. 99.98.12US).

4.2 PATIENT POSITIONING


NOTE: For the patient positioning please refer to the AMIS surgical technique (ref. 99.98.12 or ref. 99.98.12US).

4.3 SURGICAL EXPOSURE

NOTE: For the steps not described in this addendum, please refer to the AMIS surgical technique (ref. 99.98.12 or ref. 99.98.12US).


4.3.1 Skin Incision

It is crucial to understand the precise muscle position by carefully palpating the TFL (Tensor Fasciae Latae) before performing the incision. This will identify where to perform the skin incision, avoiding the Lateral Femoral Cutaneous Nerve (LFCN).



When defining the incision, the top of the greater trochanter (GT) and the Hueter's line are identified and marked.

The incision starts at the GT level and continues parallel to the skin fold up to the Hueter's line.




TIP

In line with the soft tissue sparing concept, a slightly longer skin incision should be preferred compared to a small incision with an excessive use of retractors.

The subcutaneous tissues are gently incised by the electrocautery knife along the fat lobules.

A dedicated autostatic retractor, the AMIS Bikini Beckmann, is positioned to open the incision borders.




Due to the orientation of the incision, the anterior branch of the LFCN will be visible under the skin fat layer.



Palpate the TFL again to identify the region in which it bulges the most.


The aponeurosis is incised for 1 cm and then opened with dissecting scissors in the direction of the fibers distally and proximally. The dissection runs parallel to the lateral femoral cutaneous nerve, which reduces the risk of damaging it.



9. Aponeurosis covering the rectus femoris and its incision

NOTE: With the AMIS Bikini approach, the Tensor Fasciae Latae is generally reached in a more lateral way than with the standard AMIS approach.

NOTE: Coxa vara and coxa valga hip anatomies require to shift the incision respectively proximally or distally with respect to the original incision landmark.



10. Coxa vara & Coxa valga incision shift

4.3.2 Intermuscular Approach

NOTE: For the intermuscular approach please refer to the AMIS surgical technique (ref. 99.98.12 or ref. 99.98.12US).

4.3.3 Articular Approach

NOTE: For the articular approach please refer to the AMIS surgical technique (ref. 99.98.12 or ref. 99.98.12US).

5. FEMORAL NECK OSTEOTOMY

NOTE: For the femoral neck osteotomy please refer to the AMIS surgical technique (ref. 99.98.12 or ref. 99.98.12US).

6. ACETABULAR STAGE

NOTE: For the steps not described in this addendum, please refer to the AMIS surgical technique (ref. 99.98.12 or ref. 99.98.12US).

The use of the dedicated AMIS Bikini Charnley Blades instead of the standard ones helps avoiding soft tissue injury, especially at epidermis and dermis level.




11. Dedicated AMIS Bikini Charnley Blades



12. Capsule opening with AMIS Bikini Charnley Blades

7. FEMORAL STAGE



NOTE: For the steps not described in this addendum please refer to the AMIS surgical technique (ref. 99.98.12 or ref. 99.98.12US).



The dedicated AMIS Bikini Broach Handle design has been optimized to further allow for:

- Less impingement
- Better rotational stability (control of anteversion while broaching)
- Less capsular release necessary
- Less hyperextension necessary

This results in a possible reduction of the risk of femoral fracture related to the femoral exposure.



8. REDUCTION

NOTE: For the reduction procedure please refer to the AMIS surgical technique (ref. 99.98.12 or ref. 99.98.12US).

9. WOUND CLOSURE

NOTE: For the wound closure please refer to the AMIS surgical technique (ref. 99.98.12 or ref. 99.98.12US).

10. REFERENCES

1. Menzies-Wilson, Richard & Mahalingham, Karupiah & I, Tamimi & Field, Richard. (2019) "Retrospective cohort study comparing the functional outcomes of direct anterior approach hip arthroplasty. Oblique "bikini" vs longitudinal skin incision".
2. Menzies-Wilson, Richard & Mahalingham, Karupiah & I, Tamimi & Field, Richard. (2019)." Functional Outcomes of direct anterior approach hip arthroplasty: Oblique 'bikini' versus longitudinal skin incision. 10.1177/2210491719890883.
3. Leunig, Hutmacher, Ricchiardi, Impellizzeri, Rüdiger, Naal. (2018)" Skin crease 'bikini' incision for the direct anterior approach in total hip arthroplasty: a two- to four-year comparative study in 964 patients. Bone Joint J.
4. Manrique, MD, Paskey, BS a, Tarabichi, MD, Restrepo, MD, Foltz, PhD Hozack, MD. (2019) "Total Hip Arthroplasty Through the Direct Anterior Approach Using a Bikini Incision Can Be Safely Performed in Obese Patients". J Arthroplasty

Part numbers subject to change.

NOTE FOR STERILIZATION

The instrumentation is not sterile upon delivery. Instruments must be cleaned before use and sterilized in an autoclave in accordance with the regulations of the country, EU directives where applicable, and following the instructions for use of the autoclave manufacturer. For detailed instructions, please refer to the document "Recommendations for cleaning, decontamination and sterilization of Medacta International orthopaedic devices" available at www.medacta.com.



REDEFINING BETTER
IN ORTHOPAEDICS
AND SPINE SURGERY

MEDACTA.COM



Medacta International SA
Strada Regina 34 - 6874 Castel San Pietro - Switzerland
Phone +41 91 696 60 60 - Fax +41 91 696 60 66
info@medacta.ch

Find your local dealer at: medacta.com/locations

All trademarks are property of their respective owners and are registered at least in Switzerland.
Please verify approval of the devices described in this document with your local Medacta representative.

AMIS® BIKINI
Surgical Technique Addendum

ref: 99.98.12BIKINI
rev. 01

Last update: February 2022